


## Catalog Number

9502S-12Y

**Caps, Screw Top**, with Silicone Rubber / PTFE, low durometer Septa, bonded in "knurled", polypropylene, Yellow Caps. For use with 9-425 thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

Low durometer makes the Septa easier to penetrate by the autosampler needles.

Click [HERE](#) for case quantities.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$47.55**

Printed On 2024/05/20