


Catalog Number

9502S-14C-R

Caps, Screw Top, with ultra-clean Silicone Rubber / PTFE Septa, fitted in "knurled", polypropylene, Red Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK



Additional Info

Silicone Rubber in this Septa is ultra clean and bright White.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$48.56

Printed On 2024/05/20