

## Catalog Number

9502S-50C-N

**Caps, Screw Top**, with PTFE Septa, fitted in "knurled", polypropylene, Translucent Caps. For use with 9-425mm thread, 12x32mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

Rubber-free Cap is useful when rubber cannot be used with your analysis.

This PTFE (only) Septa is for single injections or non-volatile samples as it will not reseal

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$45.84**

Printed On 2024/05/20