


Catalog Number

9502S-E-3XP

Caps, Screw Top, with Soft-Guard, White Silicone Rubber / Red PTFE Pre-Slit (X Shape) Septa, fitted in "knurled", polypropylene, Purple Caps. For use with 9-425mm thread, 12x32mm autosampler vials. Soft-Guard



Additional Info

Screw Thread: 9-425mm

Material: polypropylene, cap, silicone/PTFE septa

Features: Cross pre-slit, bonded septa

Soft-Guard cap for easier penetration of autosampler needles and resealing.

Cross slit allows needles to hit the strike target area.

Bonded septa prevents septa from falling out of the cap or into the vial

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$46.71

Printed On 2024/05/20