

## Catalog Number

9503S-05FB

**Inserts, Low Volume** for Vials, Glass. Clear, 300ul fill volume with a flat bottom and 6x31mm outer dimensions. For use in 2ml autosampler vials. BASIK Brand. (10 packs of 100). 1000/CS.



### Additional Info

First hydrolytic, Class A, Borosilicate Glass. Features a flat bottom. Residual volume is 6ul with most autosamplers.

Use with 9mm Screw Top or 11mm Snap Top 2ml, 12x32mm autosampler vials.

*Lowest cost for budget minded laboratories.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$379.06**

Printed On 2024/05/20