

Catalog Number

9503S-09-KR

Caps, Screw Top, Solid, with a Red-Orange Silicone Rubber / TEF liner, in "knurled", closed Black Caps. For use with 10-425mm thread, 12x32mm autosampler vials. MicroSolv Brand 100/PK.



Additional Info

Use with wide mouth autosampler vials. The Silicone Rubber / TEF liner provides a very good seal for longer term storage in autosampler vials.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$50.39

Printed On 2024/05/31