

## Catalog Number

9503S-0CV

**Vials, Screw Top, Glass. Clear, 2ml, 10-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Deactivated / Silanized Glass. Features a wide mouth opening for easy filling and ND10 threads.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$28.75**

Printed On 2024/05/10