

## Catalog Number

9503S-WAV

**Vials, Screw Top, Glass. Amber, 2ml, with Write-On Patch and fill lines, 10-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass, amberized. Features a wide mouth opening for easy filling and ND10 threads. Write-On Patch and fill lines are epoxy paint for reduced smearing.

Amber vials are often used for light sensitive compounds and considered low actinic.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$34.70**

Printed On 2024/05/10