


## Catalog Number

95040-40-T-D

**Caps, Screw Top**, with Silicone Rubber / PTFE, ultra-clean Septa, fitted in polypropylene, Blue Caps with a Gold, Magnetic Cover. For use with 9-425mm thread, 12x32mm autosampler vials. AQ Brand. 100/PK.



### Additional Info

Use screw top, magnetic, covered caps and avoid the difficulty in crimping and de-crimping, standard, hard metal, crimp caps. Vials with improperly crimped caps can cause evaporation and data errors; to avoid these caps that are easy to apply correctly, every time.

Pre-assembled with a magnetic gold cover over a Blue, 9-425mm, screw cap with a 1mm thick, non-slit, Silicone / PTFE Septa.

These screw caps are designed for autosamplers that use a magnetic closure such as the or the Thermo TriPlus™ and can be used with any 9mm, BASIK™, MicroSolv™ AQ™ or RSA™ Brand, Screw Top, 12x32mm, autosampler Vial.

Excellent for GC or Headspace analysis.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$118.48**

Printed On 2024/05/31