


Catalog Number

95040-SB1

Storage Box for 20ml EPA Style Vials. Violet, chemical resistant polypropylene holds 16 vials. Complete with translucent plastic lid. 1 EA.



Additional Info

Excellent for standing position on the lab bench or during transport and ideal for space saving storage in places such as refrigerators.

- Lid prevents condensation on the caps minimizing contamination during storage.
- Usable temperature range -80C to +100C.
- Box outer dimensions are 130mm wide by 130mm deep by 102mm high.
- Box and lid are Autoclavable.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$68.95

Printed On 2024/05/20