

## Catalog Number

9504S-01ND

**Inserts, Low Volume** for Vials, Glass. Clear, Deactivated, 250ul fill volume with a conical, precision point and 6x38mm outer dimensions. For use in 4ml autosampler vials. AQ Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Type 33, Deactivated / Silanized Borosilicate Glass. Features a mandrel manufactured, precision point conical tip for maximum recovery of sample and internal consistency from insert to insert. Residual volume is 3ul with most autosamplers with a 150ul useable volume. This product is manufactured in an ISO 9001:2015 Certified facility.

Springs may be required and are sold separately. Use with 13mm Screw Top 15x45mm autosampler vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$226.74**

Printed On 2024/05/20