

## Catalog Number

9504S-04ND

**Inserts, Low Volume** for Vials, Glass. Clear, Deactivated, 300ul fill volume with a conical point and 6x38mm outer dimensions with plastic spring attached. For use in 4ml autosampler vials. AQ Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Type 33, Deactivated / Silanized Borosilicate Glass. Features a conical tip for excellent performance. Residual volume is 3ul with most autosamplers with a 200ul useable volume. Polypropylene springs are attached to keep the inserts center aligned and reduce chances of the needle breaking the inserts during use. This product is manufactured in an ISO 9001:2015 Certified facility.

Use with 4ml, 13mm, Screw Top, 15x45mm autosampler vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$176.30**

Printed On 2024/05/20