

## Catalog Number

9504S-20C-BK

**Caps, Screw Top**, with PTFE / Silicone Rubber / PTFE Septa, fitted in "knurled", polypropylene, Black Caps. For use with 13-425mm thread, 15x45mm autosampler vials. MicroSolv Brand. 100/PK.



### Additional Info

Use with 15x45, 4ml autosampler vials. Used with older Waters and Shimadzu autosamplers. Designed to provide extra cleaning when "wiping" a needle, also provides a small degree of better resealing after injection compared to Silicone Rubber / PTFE Septa.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$0.00**

Printed On 2024/05/20