

## Catalog Number

9504S-35C-BLK

**Caps, Screw Top**, with Silicone Rubber / PTFE Pre-Slit Septa, fitted in "knurled", polypropylene, Black Caps. For use with 13-425mm thread, 15x45mm autosampler vials. BASIK Brand. 100/PK.



### Additional Info

Economical Pre-Slit Cap for easy penetration by autosampler needles. Recommended for older Waters® instruments. Use with 15x45mm, 4ml vials.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$37.26**

Printed On 2024/05/31