

Catalog Number

9504S-E-0CV

Vials, Screw Top, Glass. Clear, 4ml, 13-425mm threads, 15x45mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass. Features ND13 threads.

Low cost for budget minded laboratories.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$31.59

Printed On 2024/05/14