

Catalog Number

9504S-WAV

Vials, Screw Top, Glass. Amber, 4ml, with Write-On Patch and fill lines, 13-425mm threads, 15x45mm outer dimensions. For use as an autosampler vial or storage bottle. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A pharmaceutical grade, Borosilicate Glass, amberized. Features a wide mouth opening for easy filling and ND13 threads. Write-On Patch and fill lines are epoxy paint for reduced smearing.

Use these vials with Waters WISP® or older Shimadzu® autosamplers.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$53.49

Printed On 2024/05/14