

Catalog Number

9508C-10T-S

Caps, Crimp Top, with White Silicone Rubber / Red PTFE Septa, Silver Caps, with a 4mm center hole. For use with 8mm Vials. MicroSolv Brand. 100/PK.



Additional Info

These aluminum crimp caps are 8mm in diameter with a 4mm center hole and have a non-slit, Ultra-Clean, fitted, silicone rubber and PTFE septum that is 1.3mm thick. Use these caps with 8mm Serum Finish, Crimp Top Autosampler Vials.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$36.46

Printed On 2024/05/20