


Catalog Number

9508C-40T-S

Caps, Crimp Top, with 1.0mm thick, Red-Orange Natural Rubber / Transparent TEF Septa, fitted in aluminum, Silver Caps with a 4mm center hole. For use with 8mm Crimp Top autosampler vials. MicroSolv



Additional Info

This is our most popular Crimp Cap for 8mm Crimp Top Vials.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$23.23

Printed On 2024/05/20