

## Catalog Number

9508C-CB04-AV

**Vials, Crimp Top, Glass. Amber, 0.4ml. An 8mm Serum Finish with 7x30mm outer dimensions and a conical bottom. For use as an autosampler vial. MicroSolv Brand. (10 packs of 100/PK). 1000/CS**



### Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass, amberized. Amber vials are often used for light sensitive compounds and considered low actinic.

Use with 8mm aluminum crimp caps, 8mm or 9mm push caps. These conical vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

*These MicroVials have a conical outer shape and often need a holder / adapter to be used in an autosampler and will not stand by themselves.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$887.78**

Printed On 2024/05/19