

## Catalog Number

9508R-01

**Vial Rack**, for 8mm Specialty Vials. Acrylic, clear, with 50 cavities. 1 EA.



### Additional Info

For easy handling and transportation of 8mm wide autosampler vials. The clear materials makes filling easier with greater visibility compared to opaque racks. Also, racks are stackable due to the silicone feet on the base. More stable than other racks on the market due to the solid, one piece construction.

Conical and flat bottom vials can be accommodated in this 173mm wide x 95mm deep and 20mm high rack. Wells are 8.5/3mm in diameter.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$156.46**

Printed On 2024/05/20