

## Catalog Number

95096-08-CAP-M

**Caps, PE Plug Style, Polyethylene, Natural Transparent Caps.** For use with 700ul and 1.0mL, 96 well autosampler plate inserts. AQ Brand. 1000/CS.



### Additional Info

Click [HERE](#) for page with compatible inserts. Use 96 Well Plate Insert FILTER.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$113.45**

Printed On 2024/05/19