


Catalog Number

9509C-WAV-AQ

Vials, Snap Top, Glass. Amber, 2ml, with Write-On Patch and fill lines.
An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. AQ Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type-N51A Borosilicate Glass, amberized.

Compares to HPLC "Certified" vials for performance but lower priced and a great value. Bottom depths and heights are designed, and QC controlled for maximum recovery and needle safety.

Amber vials are often used for light sensitive compounds and considered low actinic.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00

Printed On 2024/05/21