


Catalog Number

9509S-10B

Caps, Screw Top, with White Silicone Rubber / White PTFE, ultra-low bleed Septa, bonded in "knurled", polypropylene, Dark Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQ Brand.



Additional Info

This very high purity Septa has low extractables and is excellent choice for HPLC. Compares to "Certified" Brands.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00

Printed On 2024/05/20