


Catalog Number

9509S-30B

Caps, Screw Top, with White Silicone Rubber / White PTFE, ultra-low bleed, Pre-Slit Septa, bonded in "knurled", polypropylene, Dark Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQ



Additional Info

Pre-slit for easy penetration by autosampler needles and *recommended for Waters® Alliance, Acquity™ and all other autosamplers that require Pre-Slit Septa.*

Good choice for LCMS & LCMS/MS.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00

Printed On 2024/05/20