

Catalog Number

9509S-WCV-AQ

Vials, Screw Top, Glass. Clear, 2ml, with Write-On Patch and fill lines, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. AQ Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type-33 Borosilicate Glass. Write-On Patch and fill lines are ceramic and will resist smearing or scratching when working with solvents.

Compares to HPLC "Certified" vials for performance but lower priced and a great value. Bottom depths and heights are designed, and QC controlled for maximum recovery and needle safety.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$37.85

Printed On 2024/05/10