

## Catalog Number

9512C-0CV-P

**Vials, Snap Top, Glass. Clear, 1.5ml, U Shape. An 11mm snap-ring with 12x32mm outer dimensions and a fluted stand. For use as an autosampler vial. MicroSolv Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass.

Compares to Champagne Style Vials. Features a 15ul residual volume for most autosamplers due to the center draining, inner U shape.

*We recommend to use these vials with [Soft-Guard™](#) Snap Caps.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$80.01**

Printed On 2024/05/21