

Catalog Number

9512S-3MP-RS

Vial & Cap Kit Includes 100 EA 1.2ml, MRQ, Screw Top, Clear RSA™ Autosampler Vials and 100 EA matching Screw Caps with Clear AQR Silicone Rubber / Clear PTFE, ultra-pure, Pre-Slit Septa, fitted in the



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Reduced Surface Activity Vials: Cat. No [9512S-0CV-T-RS](#) and LCMS/MS compatible Caps: Cat No. [9509S-30C-B](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$257.58

Printed On 2024/05/20