

## Catalog Number

9516S-0CV

**Vials, Screw Top, Glass. Clear, 5ml, Max Recovery™, 16mm threads, 22x37mm outer dimensions. For use as a autosampler vial or storage bottle. AQ Brand. 30/PK.**



### Additional Info

First hydrolytic, Class A, Type-33 Borosilicate Glass. This market leading 5ml vial (total fill volume 6ml) has a low residual volume due to the center-draining, conical bottom, which allows needles with bottom entry ports to recover the maximum amount of liquid during use.

Bottom depths are designed and QC controlled to prevent needles from hitting the bottom of the vials. This vial with the proper cap does not exceed 40mm.

Click [HERE](#) for compatible Screw Caps.

*Compare to the high recovery vial used in the Agilent InfinityLab™ LC Series; fully compatible.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$195.70**

Printed On 2024/05/13