


## Catalog Number

95296-10P-EA

**Mat, Sealing,** Silicone Rubber, Plug Style for 96, 350ul Inserts. U-2D  
Brand. 1 EA.



### Additional Info

Made of ultra pure Silicone Rubber, this mat is designed and molded to seal up to 96 each 350ul Inserts in the U-2D Micro-Sample Management System. This mat is not Pre-Slit and has alphanumeric markings on the top side of the mat, along the edges for identification of the Inserts.

Click [HERE](#) for case quantity.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$19.73**

Printed On 2024/06/01