


## Catalog Number

95296-30P-CS

**Mats, Sealing,** Silicone Rubber, Pre-Slit, Plug Style for 96, 350ul Inserts.  
U-2D Brand. 5/CS.



### Additional Info

Made of ultra pure Silicone Rubber, this mat is designed and molded to seal up to 96 each 350ul Inserts in the U-2D Micro-Sample Management System. This mat is Pre-Slit and has alphanumeric markings on the top side of the mat, along the edges for identification of the Inserts.

Click [HERE](#) for individual pack quantities.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$125.26**

Printed On 2024/05/20