

Catalog Number

95296-350-PK

Inserts, in Racks. Includes 96, Clear, 350ul Glass Low Volume Precision Inserts with wide conical points. For use in the Micro-Sample Management System. U-2D Brand. 6/PK.



Additional Info

Inserts are 1st Hydrolytic, Class A, TYPE 33 Borosilicate Glass made for use with the U-2D Micro-Sample Management System or in 2ml, 12x32mm AQ Brand Autosampler Vials.

These Inserts have a wide conical point for fast sample filling which minimizes bubble formation. Pre-Loaded in Racks.

Click [HERE](#) for case quantity. *Bases, Sealing Film or Mats are sold separately.*

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$651.73

Printed On 2024/05/20