

## Catalog Number

95296-MARK-F5

**Marker**, Laboratory, permanent. Ultra-Fine tip, 0.45mm. Black Ink. 1 EA.



### Additional Info

These specially made markers are perfect for making permanent, Ultra-Fine line details in micro-spaces or on unusual/challenging surfaces. Excellent for labelling vials, U-2D MicroSample Management Systems and other laboratory devices.

These "pens" can be used to make precise, long-lasting marks wherever there is a need for waterproof ink that will write on non-porous surfaces.

Laboratory Markers features:

- A fine line nib with a protected tip that will not split.
- Chemically stable and waterproof alcohol-based ink.
- Low odor and fulfills non-toxicity requirements.
- Permanent on most surfaces.
- Alcohol based cleaners can remove ink from non-porous surfaces.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$10.90**

Printed On 2024/05/20