

## Catalog Number

9532C-0XV-K

**Vials, Snap Top, Plastic.** Clear, 300ul MicroVial, TPX. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 1000/CS.



### Additional Info

This narrow mouth, vial is manufactured with virgin, medical grade polymethylpentene which produces the visually clearest, translucent, built-in insert with a 300ul total fill capacity. Residual volume can be as low as 2ul depending on the autosampler and viscosity of samples used.

Popular for samples that react with glass and often used in bio-analysis with aqueous samples. Reactivity is similar to polypropylene but is more expensive and chosen when seeing the sample in the vial is important.

Click [HERE](#) for smaller pack quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$592.78**

Printed On 2024/05/21