

## Catalog Number

9532C-TS-A

**Vials, Snap Top, Plastic. Amber / Clear Glass Insert, 200ul Fused Insert™.** An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.



### Additional Info

These narrow mouth vials are made with amberized, TPX (polymethylpentene) Plastic and have a first hydrolytic, Type N51A, Clear, pharmaceutical grade, Borosilicate Glass Insert, with a conical bottom, embedded in the plastic. These vials have a 200ul total fill volume and will not leak from the bottom if the needle crashes the insert. The Insert will not move in the vial.

Amber vials are often chosen for light sensitive compounds and are considered low actinic. Compares to TopSert, GlassPlastic™ and Glastic™ Brands.

Click [HERE](#) for Clear Vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$237.73**

Printed On 2024/05/17