

## Catalog Number

9532S-06-0DCP

**Vial & Cap Kit** Includes 100 EA 300ul, Fused Insert, Screw Top, Clear Deactivated Autosampler Vials with Write on Patch and fill lines and 100 EA matching Screw Caps with Clear Silicone Rubber / Clear Tefzel, ultra-



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Silanized / Deactivated Vials: Cat. No [9532S-0DCV](#) and Caps: Cat No. [9502S-06-LCC](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$277.85**

Printed On 2024/05/20