

## Catalog Number

9532S-1CP-WT-RS

**Vial & Cap Kit** Includes 100 EA 300ul, Fused Insert with a wide tip, Screw Top, Clear RSA™ Autosampler Vials with Write on Patch and fill lines and 100 EA matching Screw Caps with Clear AQR Silicone Rubber



### Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Reduced Surface Activity Vials: Cat. No. [9532S-WCV-WT-RS](#) and Caps: Cat No. [9509S-10C-B](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$250.63**

Printed On 2024/05/19