

Catalog Number

9532S-3AP-RS

Vial & Cap Kit Includes 100 EA 300ul, Fused Insert, Screw Top, Amber RSA Autosampler Vials with Write on Patch and fill lines and 100 EA matching Screw Caps with Clear AQR Silicone Rubber / Clear PTFE,



Additional Info

Easy Purchase Packs are a convenient way to order and inventory autosampler vials with their matching caps. The kit packaging will open and re-close making it useful for keeping vial and cap types matched, organized and clean during use.

Click on links below for more Specifications for the 9mm Vials and Caps:

Supplied with Reduced Surface Activity Vials: Cat. No. [9532S-0AV-RS](#) and Caps: Cat No. [9509S-30C-B](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$233.04

Printed On 2024/05/19