

## Catalog Number

9532S-TS-C

**Vials, Screw Top, Plastic.** Clear / Clear Glass Insert, 200ul Fused Insert™, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage container. MicroSolv Brand. 100/PK.



### Additional Info

These narrow mouth vials are made with TPX (polymethylpentene) Plastic, and have a first hydrolytic, Type N51A, Clear, pharmaceutical grade, Borosilicate Glass Insert, with a conical bottom embedded in the plastic. These vials have a 200ul total fill volume and will not leak from the bottom if the needle crashes the glass insert. The insert will not move in the vial.

*Compares to TopSert, GlassPlastic™ and Glastic™ Brands.*

Click [HERE](#) for Amber Vials.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$228.34**

Printed On 2024/05/12