

## Catalog Number

9534S-E-0CV-K

**Vials, Screw Top, Glass. Clear, 3.8ml, Max Recovery™, 13-425mm threads, 14.75x45mm outer dimensions for 13mm screw caps. For use as an autosampler vial or storage container. BASIK Brand. 1000/CS.**



### Additional Info

First hydrolytic, Class A, pharmaceutical grade, Borosilicate Glass. This 3.8ml vial has the low residual volume of 18ul.

*Low priced for budget minded laboratories.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,703.99**

Printed On 2024/05/13