

## Catalog Number

9572S-MS

**Vials, Screw Top, Plastic.** Clear, 700ul MicroVial, LCMS Compatible, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. AQ Brand. 100/PK.



### Additional Info

This narrow mouth vial is manufactured with ultra-pure, virgin, medical grade polypropylene and is translucent with a 700ul total fill capacity. Residual volume can be as low as 75ul depending on the autosampler and viscosity of samples used.

This specially formulated plastic is used in our special molding process where nothing is added to the process producing a very unique, plastic vial is that compatible with LCMC (extractables free in most LCMS solvents). This product is manufactured in an ISO 9001:2015 Certified facility.

Popular for samples that react with glass and often used in bio-analysis with aqueous samples.

*Excellent for biological compounds with LCMS.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$65.69**

Printed On 2024/05/09