

## Catalog Number

9572S-PP-K

**Vials, Screw Top, Plastic.** 700ul, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 1000/CS.



### Additional Info

These injection molded, translucent, plastic vials are made with virgin polypropylene and chosen for bio-analytical reasons.

Click [HERE](#) for smaller pack quantities.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$302.60**

Printed On 2024/05/13