

Catalog Number

9582S-03-B-M

Caps, Screw Top, with Silicone Rubber / PTFE Pre-Slit Septa, fitted in "knurled", polypropylene, Black Caps. For use with 8-425mm thread, 12x32mm autosampler vials. BASIK Brand. 1000/CS.



Additional Info

Wider hole for use with narrow mouth, 2ml vials. Wider hole provides a larger target area for the autosampler needle during injection.

Click [HERE](#) for smaller pack quantities

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$317.61

Printed On 2024/05/20