

## The Vial Centrifuge

CE10-C0817-VE

Part No.

**The Vial Centrifuge** with 8 Position rotor for 2ml heavy wall vials, Includes 8 special adapters one pack of RSA-Pro X™ MRQ screw top vials with black, non-slit AQR™ Caps. (9512S-1DMP-PX). MicroSolv™



### Additional Info

High Speed Centrifuge with a specially designed rotor and 8 adapters and cushions for **approved** glass MRQ™ (1.2ml) and Max Recovery (1.8ml) single piece, heavy wall autosampler vials. Includes one pack of [RSA-Pro X™](#) Hydrophobic, MRQ™ screw top vials with black, non-slit AQR™ high purity screw caps (9512S-1DMP-PX), and a wrench for tightening the rotor in the centrifuge.

The Vial Centrifuge offers 16,800xg with a unique air-flow system that keeps your samples cool. This unit features a large, touch panel screen for easy operation which controls and displays actual speed (in rpm or rcf) along with set and remaining spin times. Store up to 6 speed/time profiles in memory as the home screen displays a "mySpin" icon for instant recall. The mySpin™ button program recall makes it easy to get results fast.

The Vial Centrifuge is cold room safe.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$3,810.24**

Printed On 2024/05/20